

# User's manual



# *I-Press 6C Pressure Therapy Device*

## Instructions for use & Technical description

**Please read these instructions carefully before using  
your new device!  
This manual is an integral part of the equipment and  
must be kept until it is destroyed.**

**This equipment has been designed and manufactured  
for therapeutic use.  
The use is reserved to professionals Physiotherapists**

**If you have a problem or do not understand this manual,  
please contact your distributor (see stamp on the last page) or  
contact Électronique du Mazet at :**

**Tel: (33) 4 71 65 02 16 - Fax: (33) 4 71 65 06 55**

**Please return the warranty certificate on the last page of this  
manual within 15 days of installation or receipt.**

# 1 Presentation of the device

The i-Press 6C is a pressure therapy device that can be used in physiotherapy to help treat venous disorders and lymphoedema. It can also be used for wellness programmes.

The computerised technology used in the i-Press 6C allows for ease of use and easy menu navigation.

The programs, predefined in the device, allow most of the pneumatic drainage techniques to be performed:

- 1. venous system: venous insufficiency, varicose veins, water retention problems "heavy legs", ulcers, prevention of deep vein thrombosis (DVT).
- 2. lymphatic system :
  - Lower Extremity** Zone Lymphatic System Dysfunction: Secondary Lymphoedema, Edema
  - Upper Extremity** Zone: Lymphedema
- Well-being: Relaxation (non-therapeutic programme will not be covered in the clinical evaluation), Recovery from Exercise

For all these programmes, the parameters that can be changed are as follows:

- Exerted pressure (mmHg)
- Duration of treatment
- Treated area
- Working time
- Rest time

These settings can be changed and saved.

The i-Press 6C offers a choice of 2 operating modes:

- Access to treatment through a clinical guide according to pathologies with preset, but modifiable and recordable parameters  
This option allows for ease of use and provides a guarantee of safety for the user.
- Access through the customised processing base where all parameters can be modified and saved.  
This option allows the programme to be adapted to a particular need.

## Table of Contents :

1	Introduction to the device.....	3
2	Description and technical information .....	5
2.1	Symbols used.....	6
2.2	Technical characteristics.....	7
2.2.1	General characteristics.....	7
2.2.2	Technical characteristics of the device .....	7
2.2.3	Accessories.....	8
2.2.1	Applied parts .....	8
2.3	Nameplate label .....	9
2.3.1	Device type label .....	9
2.3.2	Accessory type label.....	9
2.4	Warnings.....	10
2.5	Precautions.....	11
2.5.1	Residual risks .....	11
2.5.2	Mains failure .....	11
2.5.3	Treated areas .....	11
2.6	Confidentiality of patient data .....	11
3	Installation of the appliance .....	12
4	User's Manual.....	13
4.1	Getting started with the device .....	13
4.1.1	Powering up / starting / stopping .....	13
4.1.2	Using the touch screen.....	13
4.1.3	Main menu .....	14
4.2	Choosing a treatment .....	15
4.2.1	From the last treatment performed.....	15
4.2.2	From the diagnosis .....	15
4.2.3	From customised programmes .....	15
4.3	Changing the parameters .....	16
4.4	Performing a treatment .....	16
4.4.1	Starting the treatment .....	16
4.4.2	During treatment .....	17
4.4.3	End of treatment.....	18
4.5	Saving a treatment .....	18
4.6	Technical Information, Configuration and Settings .....	19
5	Clinical Guide .....	20
5.1	Target population.....	20
5.2	Expected performance .....	20
5.3	Major contraindications .....	21
5.4	Side-effects .....	21
6	Maintenance, servicing.....	22
6.1	Enclosure .....	22
6.2	Accessories .....	22
6.3	Sterilisation :.....	22
7	Malfunction .....	23
8	After-sales service and warranty .....	25
9	Disposal.....	26
9.1	Accessories .....	26
9.2	Electronics .....	26
10	Transport and storage .....	26
11	EC Declaration .....	27
12	Manufacturer .....	27
13	EMC compliance table .....	28

## 2 Description and technical information

-This Operation and Maintenance Manual is published to help you get to grips with your i-Press 6C from the initial acceptance and commissioning stage through to subsequent operation and maintenance.

If you have any difficulty in understanding this manual, contact the manufacturer, Électronique du Mazet, your dealer or distributor.

-This document must be kept in a safe place, protected from atmospheric agents, where it cannot be damaged.

This document guarantees that the devices and their documentation are technically up-to-date at the time of marketing. However, we reserve the right to make changes to the device and its documentation without any obligation to update these documents.

-In the case of transfer of the device to a third party, it is obligatory to inform Électronique du Mazet of the details of the new owner of the device. It is imperative to provide the new owner with all documents, accessories and packaging relating to the appliance.

-Only personnel who have been informed of the contents of this document may use the equipment. Failure to comply with any of the instructions contained in this document will release Électronique du Mazet and its authorised distributors from the consequences of accidents or damage to personnel or third parties (including patients).

## 2.1 Symbols used



**Warning:** this logo draws your attention to a specific point



**Operating instructions:** this logo informs you that the operating instructions must be read for safe use of the appliance




**Type B applied part:** applied part in contact with the patient and which can be connected to earth.



**Recycling:** This appliance should be disposed of at an appropriate collection and recycling facility. Consult the manufacturer.



**Protective earth**



**Fuse**



**Caution: Switching the** appliance off and on



Alternating current



Serial number



Manufacturer



Date of manufacture



Country of manufacture



Product reference

## ***2.2 Technical specifications***

### **2.2.1 General characteristics**

- Operating temperature: 0°C to 40°C.
- Storage temperature: -40°C to 70°C.
- Operating relative humidity: 30% to 75%.
- Operating altitude: < 2000 metres
- Operating pressure: between 80 and 110 kPa

### **2.2.2 Technical characteristics of the device**

- Housing dimensions: 340 x 320 x 140 mm
- Case weight: 3.3 Kg
- Housing colour: white - metallic grey screen
  
- Power supply: 230VAC - 50Hz
- Power consumption: <50VA (230VAC)
- Fuses: 2x size 5x20mm - T2AH-250V
- Class I electrical equipment
- Power on indication: Display illumination
  
- Medical Class IIa equipment.
- Type B applied part

### 2.2.3 Accessories

This device is delivered with the following accessories as standard:

- 1 Power cord
- 1 Pair of 6-cell boots
- 1 Blanking plug
- 1 User's Manual
- 1 Clinical Guide

The optional accessories available are :

- 5-cell sleeve
- 5-cell belt

The use of accessories not recommended by the manufacturer does not engage his responsibility

### 2.2.1 Applied parts










Treatment accessories must not be placed in direct contact with the patient's skin. A single-use hygiene sleeve must be used, which is considered a **Type B applied part**. Hygiene sleeves are not supplied with the appliance.





## 2.3 Nameplate label

The information and characteristics are given on the back of each appliance on a label.

### 2.3.1 Device nameplate




<b>I-Press 6C</b>		<b>REF</b> MEG009KP101-A0
<b>SN</b> PM06-0155		<b>2020</b>
 230VAC 50VA 50Hz	 Rating T2AH-250V 2x size 5x20	
	<b>Electronique du Mazet</b> ZA Route de Tence 43520 LE MAZET St VOY	 FRANCE
		
	 Type B	<b>CE</b> 0459




	230VAC 50VA
Hz	50Hz
	Rating T2AH-250V 2x size 5x20



### 2.3.2 Accessory identification label

<b>SN</b> XXXXXX-XXXX
-----------------------

<b>5P sleeve</b>		<b>CE</b>
<b>REF</b> MEG009KP502A0		
	ELECTRONICS of MAZET ZA Route de Tence 43520 LE MAZET St VOY	 FRANCE

<b>Boot 6P</b>		<b>CE</b>
<b>REF</b> MEG009KP503A0		
	ELECTRONICS of MAZET ZA Route de Tence 43520 LE MAZET St VOY	 FRANCE

## 2.4 Warnings



**CAUTION**: Install the appliance on a flat, stable surface. Do not block any ventilation openings (no objects closer than 4cm).



**CAUTION**: Power strips must not be placed on the floor. No other electrical appliance or power strip should be connected to the power strip.



**CAUTION**: The appliance must be connected to an outlet with a grounding terminal (Class I electrical appliance)



**CAUTION**: The unit must be positioned so that the mains cable is freely accessible in case of emergency.



**CAUTION**: In case of emergency, disconnect the mains cable directly from the unit.



**CAUTION**: No modifications to the device are permitted. It is strictly forbidden to open the housing of the device.



**CAUTION**: This equipment complies with applicable electromagnetic compatibility standards. If you experience interference or other problems with another device, contact Électronique du Mazet or the distributor for advice on how to avoid or minimize the problem.



**CAUTION**: This equipment is not intended for use in residential environments and may not provide adequate protection for radio reception in such environments.



**CAUTION**: Operating altitude below 2000m.  
The performance of the aircraft decreases with altitude.



**CAUTION**: The appliance must be used with the accessories supplied by the manufacturer.



**CAUTION**: The device must not be accessible to the patient.  
It should not be placed in contact with the patient.

## ***2.5 Precautions***

### **2.5.1 Residual risks**

### **2.5.2 Mains failure**

In the event of a power failure during treatment, it is advisable to disconnect the pneumatic connectors from the back of the device to relieve pressure on the limbs.

### **2.5.3 Treated areas**

The areas treated by the device are the lower or upper limbs.  
See **§7-clinical guide** for more information.

## ***2.6 Confidentiality of patient data***

The device collects data when a treatment is saved. The data is stored in the device. It is the responsibility of the practitioner to apply and comply with the European Parliament's General Data Protection Regulation 2016/679.

When returning to the After Sales Service, the practitioner must delete the patient data so that it is not disclosed.

### 3 Installation of the appliance

Open the carton, remove the accessories and the i-Press 6C

Remove the lightweight plastic wrappings covering the unit.  
Check the contents of the box against the **packing list** included with the documentation.

Check that the contents of the box are not damaged; if you have any doubts about the integrity of the appliance or its accessories and that the correct operation of the appliance could be questioned, contact Électronique du Mazet


If the appliance was stored in a cold place and there was a risk of condensation, **let the appliance stand for at least 2 hours at room temperature** before switching on.

Before using the appliance for the first time, it is advisable to clean it and its accessories (see **§8 Maintenance**).

Install the device on a stable stand at working height and out of the patient's environment.



-Connecting the power cord  
-Toggle the switch :  
Position 0: Off  
Position 1: On



-Connect the accessories: Clip on the connectors, respecting the coding.  
-Disconnecting accessories :  
Press the blue button and pull



Example of connection of 2 accessories:  
**Boots**



Example of connection 1 accessory:  
**Sleeve**  
In this case the unused track must be closed with a plug (supplied)

## 4 User's manual

### 4.1 Handling the device

#### 4.1.1 Power on / start / stop

Turn on the power using the switch on the back of the unit (Position I: On / position 0: Off) (See §5 Installation of the appliance). The home screen lights up and displays the software version.



Choice of display language

#### 4.1.2 Use of the touch screen

The lists of choices displayed on the screen and the validations and navigation in the menus are indicated by "action" buttons on the touch screen. To access the desired function, press in the indicated area.


### 4.1.3 Menu main

By pressing the corresponding key, this menu allows access to :

- To the "**Pathology Access**" treatment database
- To personalised treatments: "**Personalised Base**"  
(See §6.2.3 Choosing a treatment from the personalised programmes)
- To technical information and settings: "**Configuration Settings**" (See §6.6)



All the menus accessible afterwards will be equipped with a button allowing the return to the main menu (placed at the top left)



## **4.2 Choice of treatment**

### **4.2.1 From the last treatment performed**


To select the last treatment performed, press the bottom left screen when the start page is displayed.

### **4.2.2 From u diagnosis**

Selecting the pathology to be treated, by pressing the corresponding action button, allows access to the menu **§6.3 Modification of parameters**.

Details of the predefined parameters for each condition are given in **§7 Clinical Guide**.

### **4.2.3 From the customised programmes**



Your own programs (up to 40 slots), distinguished by the name you assign to them (16 characters available), can be assigned to some of your patients or to specific pathologies .

Press on the selected treatment to access the menu **§6.3 Modification of parameters**.


### **4.3 Modification of parameters**

When accessing the treatment, and before starting it, it is possible to modify the "Time" (treatment) and "Pressure" (in the cells) parameters by pressing the "+" or "-" key of the one you wish to set.


### **4.4 Performing a treatment**

#### **4.4.1 Start treatment**

Choose the category of treatment you want.




Choose the desired treatment.




Once the treatment has been chosen, adapt the parameters to the patient and press confirm.





This icon gives access to the treatment settings:



-**Working time** from 0 to 15 seconds  
(keeping the cell pressurised)

-**Resting time** of 4 to 15 seconds  
(Deflation between 2 inflations)

- **Number of accessories:** (setting available with leg programs) Default set to 2.

- **Stopping during the cycle :**


-Valid: end of treatment at the end of the timer.

-Otherwise the programme stops at the end of a complete cycle.

#### 4.4.2 During treatment

All parameters can be changed during processing by selection on the touch screen.

The treatment can be stopped or restarted by selection on the touch screen.



### 4.4.3 End of treatment

End of treatment.




The end of the treatment is signalled by a succession of beeps.

You can stop the active deflation to return to the general menu

### 4.5 Saving a treatment

After stopping, at the end of a treatment or before starting the treatment, the user has the possibility to save the parameters of his treatment in one of the programmes of the personalised database (see § 6.2.3).

This is done by pressing the "save" button from the processing start screen.



Choosing a programme for the custom base :


-Enter the name of the treatment (max. 16 characters) for a total of 40 records.

## **4.6 Technical Information, Configuration and Settings**

This screen provides access to hardware technical information, menu language selection, screen brightness adjustment, sound information selection and testing.



You will then have access to the following menu:



- **The following information** gives you our contact details and those of our after-sales service
- **Language:** selects the language of the device (English, French or Spanish)
- **Screen:** allows you to adjust the contrast
- **Beep sound:** allows you to activate or deactivate the beep sound when a key is pressed (the beep at the end of the treatment cannot be deactivated)
- **Touch Panel:** allows you to adjust the sensitivity of the touch panel
- **Reserved access:** allows you to launch a self-diagnosis of the device using a code that will be communicated to you by our after-sales service (Information button) in the event of a failure or malfunction of the device.

## 5 Clinical Guide

### **5.1 Target population**

The device is intended for use by adults of any gender over the age of 18.

### **5.2 Expected performance**

The **i-Press 6C** is a pressotherapy device to assist in the management of venous and lymphatic pathologies by performing mechanical lymphatic drainage.

This technique allows to favour the so-called "return" circulation. Indeed, the interest of lymphatic drainage is no longer to be demonstrated. The lymphatic system is a major function of our body. The transport of the lymph is not always ensured properly, moreover, it is rather slow. This can have many repercussions on the circulatory system, as well as on the general state of the body. When we know the importance of the lymph in the good functioning of our defence system and the elimination of toxins, we understand better the interest of optimising lymphatic drainage. Commonly used, manual drainage (MD) is practiced by a good number of physiotherapists. This technique is very widespread and does not always offer all the keys to the various lymphatic pathologies. A specific case of chronic oedema shows the need to treat it in synergy with other techniques: bandaging and pneumatic drainage (PD).

Pneumatic drainage must come as a complement, not as a replacement of manual techniques. It brings to the practitioners an additional help to the treatment of their patients, even an essential tool in the treatment of oedemas. Pneumatic drainage leads to decongestion, it will then be used before a manual drainage, while in other cases it can be used simultaneously, or even after a manual drainage to reinforce and prolong its effects.

Thanks to the new possibilities of programming and adjustment which allow numerous movements, we speak of sequential and programmable pressotherapy, the indications for the use of pneumatic drainage are numerous and can be classified into 3 categories:

1. **Venous system:** Venous insufficiency, Varicose veins, Water retention problems "heavy legs", Ulcers, Prevention of deep vein thrombosis
2. **Lymphatic system:**  
**Lower Extremity** Zone Lymphatic System Dysfunction: Secondary Lymphoedema, Edema  
**Upper Extremity** Zone: Lymphodema
3. **Well-being:** Relaxation, Recovery after exercise

### ***5.3 Major contraindications***

This appliance **must not be used** in the following cases:

- Emboligenic deep vein thrombosis
- Untreated heart failure, Severe arterial disease
- Skin infection (Eresypelas, urticaria, ...)
- Sensitivity disorder
- Severe renal failure
- Febrile state (fever, ...)
- Lymphangitis
- Systemic oedema

**Note:** Pressotherapy treatments for pregnant or breastfeeding women have not been clinically studied, therefore treatment of these individuals is the responsibility of the practitioner.

The contraindications are not exhaustive and we advise the user to seek advice in case of doubt.

### ***5.4 Side effects***

To date, the medical literature does not mention any side effects of pressotherapy.

## 6 Maintenance, servicing

The I-press 6C is designed to have a life of 5 years.

**To ensure that the performance of the device is maintained throughout its life, it is necessary to have the device checked by Électronique du Mazet technicians every 2 years.**

### ***6.1 Housing***

The casing only requires normal, periodic cleaning of its outer surface, which may become dirty. The same applies to the power cord.

The touch screen should be cleaned with a soft, dry cloth, **without any product or water.**

Clean the rest of the appliance only with a dry or slightly damp cloth.  
Be sure to unplug the power cord before cleaning.

### ***6.2 Accessories***

**The treatment accessories must not be placed in direct contact with the patient's skin. A hygiene sheath must be used.**

The treatment accessories can be cleaned with a dry or slightly damp cloth.

### ***6.3 Sterilization :***

This device is not sterile,  
The accessories are not sterile, nor are they intended to be sterilised.

## 7 Malfunction

If you notice a malfunction that is not commented on in the documents accompanying the appliance (see below), please inform your distributor or the manufacturer.

In case of a shipment of the device, please observe the following instructions:

- Decontaminate and clean the unit and its accessories.
- Use the original packaging, including the retaining flanges.
- Attach all accessories to the appliance.
- Set up the various elements.
- Ensure that the packaging is properly sealed.

Shipping address :

**Electronique du Mazet  
ZA Route de Tence  
43520 Le Mazet St Voy**

**Tel: (33) 4 71 65 02 16**

**Fax: (33) 4 71 65 06 55**

**E-mail: [sav@electroniquedumazet.com](mailto:sav@electroniquedumazet.com)**

Possible malfunctions:

<b>Description of the anomaly</b>	<b>Possible causes</b>	<b>Actions</b>
Screen off	Problem with the electricity network	Check the mains connection
	Switching on the appliance	Check the position of the on/off switch (position I)
	Fuses out of order	Check and change fuses
	Other cause	Contact the After Sales Service
Pockets do not bulge	Poorly connected accessory	Check the locking of the pneumatic connectors
	Defective accessory (leakage)	Contact the After Sales Service
	Pneumatic problem (faulty pump or solenoid valve)	Contact the After Sales Service
Pockets do not deflate (defective active deflation)	Pneumatic problem (faulty pump or solenoid valve)	Contact the After Sales Service
Other anomaly	Unknown	Contact the After Sales Service

If the device is dropped or if water penetrates, it is imperative to have the device checked by Électronique du Mazet to exclude any risk (patient and user) related to the use of the device.



## 8 After-sales service and warranty

This appliance is guaranteed by your supplier under the conditions specified in this document, provided that :

- Only accessories supplied by Électronique du Mazet or its distributors should be used.
- Any modification, repair, extension, adaptation and adjustment of the appliance must be carried out by Électronique du Mazet or its authorised distributors for these operations.
- The working environment meets all regulatory and legal requirements.
- The appliance may only be used by competent and qualified personnel. Use must be in accordance with the instructions in this user's manual.
- Treatments should only be used for the applications for which they are intended and which are described in this manual.
- The appliance must be regularly maintained according to the manufacturer's instructions.
- All legal requirements for the use of this device are met.
- The appliance uses only the accessories supplied or specified by the manufacturer.
- Machine parts and spare parts must not be replaced by the user.

Inappropriate use of this device or neglect of maintenance relieves Électronique du Mazet and its authorised distributors of all responsibility for defects, breakdowns, malfunctions, damage, injuries and the like.

The warranty is void if the operating instructions in this manual are not strictly followed.

**The warranty period is 24 months from the date of delivery of the device.**

**The accessories are guaranteed for 6 months from the date of delivery of the device.**

**Transport and packaging costs are not included in the guarantee.**

## **9 Disposal**

### ***9.1 Accessories***

As soon as any damage to an accessory is detected, the product must be cleaned with a broad spectrum disinfectant and returned to the manufacturer.

### ***9.2 Electronics***

Should the I-Press 6C device fail to function or become unusable, please return it to the manufacturer or take it to a Recylum collection point.

As part of its commitment to the environment, Électronique du Mazet finances the Réylum recycling network dedicated to WEEE Pro, which takes back free of charge electrical lighting equipment, control and monitoring equipment and used medical devices (more information on [www.recylum.com](http://www.recylum.com)).



## **10 Transport and storage**

The appliance must be transported and stored in its original packaging or in packaging that protects it from external damage.

Store in a clean, dry place at room temperature.

## 11 CE declaration

ÉLECTRONIQUE DU MAZET can provide the CE declaration for this device on request.

The first affixing of the medical CE on this device took place on 14/10/2016

## 12 Manufacturer

Électronique du Mazet is a company located in the heart of the Massif Central. Originally a simple manufacturer of electronic cards, over the years it has developed its own brand of medical equipment, mainly for physiotherapy.

Today, EDM studies, develops, manufactures and markets pressotherapy, depressotherapy and electrotherapy (urological rehabilitation) equipment.

For further information, please do not hesitate to contact us.

***SAS Électronique du Mazet***  
***ZA Route de Tence***  
***43520 Le Mazet St Voy***  
**Tel: +33 (0)4 71 65 02 16**  
**Fax: +33 (0)4 71 65 06 55**




[facebook.com/mazet-santé](https://facebook.com/mazet-santé)



[www.electroniquedumazet.com](http://www.electroniquedumazet.com)

## 13 EMC compliance table

<b>EMC compliance to IEC / EN 60601-1-2 (2014)</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or user of the equipment should ensure that it is used in such an environment.			
<b>Emissions testing</b>	<b>Standard</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidelines</b>
RF emissions	CISPR 11	Group 1	The device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause interference in any nearby electronic device.
RF emissions	CISPR 11	Class B	The appliance is suitable for use in all premises, including domestic premises and those directly connected to the public low-voltage power supply network supplying buildings for domestic use.
Harmonic emissions	IEC 61000-3-2		Class A
Voltage fluctuations / Flicker	IEC 61000-3-3		Compliant
The device is intended for use in the electromagnetic environment specified below. The customer or user of the equipment should ensure that it is used in such an environment.			
<b>IMMUNITY test</b>	<b>Test level IEC 60601</b>	<b>Level of compliance</b>	<b>Electromagnetic environment - guidelines</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic materials, the relative humidity should be at least 30%.
Transients fast in bursts IEC 61000-4-4	± 2 kV for lines power supply electric ± 1 kV for lines input/output	± 2 kV for lines power supply electric	The quality of the power supply should be that of a typical commercial or hospital environment.
Surge voltage transitional IEC 61000-4-5	± 1 kV between phases ± 2 kV between phase and earth	± 1 kV between phases ± 2 kV between phase and earth	The quality of the power supply should be that of a typical commercial or hospital environment.
Tension dips, short cuts and variations in tension on input lines power supply electric IEC 61000-4-11	<5% TU (>95% UT trough) for 0.5 cycle 40% UT (60% of UT troughs) for 5 cycles 70% UT (30% of UT trough) for 25 cycles <5% TU (>95% UT trough) for 5 s	<5% TU (>95% UT trough) for 0.5 cycle 40% UT (60% of UT trough) for 5 cycles 70% UT (30% of UT trough) for 25 cycles <5% TU (>95% UT trough) for 5 s	The quality of the power supply should be that of a typical commercial or hospital environment. If the user of the device requires continuous operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery. NOTE UT is the AC mains voltage before the test level is applied.
Magnetic field at the frequency of the electrical network (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the frequency of the power system should have levels characteristic of a representative location in a typical commercial or hospital environment.
<b>IMMUNITY test</b>	<b>Test level IEC 60601</b>	<b>Level of compliance</b>	<b>Electromagnetic environment - guidelines</b>

<p>RF disturbances conducted IEC 61000-4-6</p>	<p>3 Vrms 150kHz-80MHz 3V/m 80MHz-2.5GHz</p>	<p>3 Vrms 3V/m</p>	<p>Portable and mobile RF communications equipment should not be used closer to any part of the equipment, including cables, than the recommended separation distance, calculated from the equation for the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.67 \cdot \sqrt{P}$ $d = 1.67 \cdot \sqrt{P} \text{ } 80\text{MHz}-800\text{MHz}$ $d = 2.33 \cdot \sqrt{P} \text{ } 800\text{MHz}-2.5\text{GHz}$ <p>where <math>P</math> is the maximum output power characteristic of the transmitter in watts (W), according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>The field strengths of fixed RF transmitters, as determined by an on-site electromagnetic investigation, should be below the compliance level in each frequencies. b</p> <p>Interference may occur in the vicinity of the device marked with the following symbol:</p> 
--	--	------------------------	---

NOTE 1 At 80 MHz and 800 MHz, the highest frequency range applies.  
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflections from structures, objects and people.

- a) The field strengths of fixed transmitters, such as base stations for radiotelephones (cellular/wireless) and land mobile radios, amateur radio, AM and FM broadcasting, and TV broadcasting, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an on-site electromagnetic investigation should be considered. If the field strength, measured at the location where the equipment is used, exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be required, such as reorienting or repositioning the equipment.
- b) In the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the device**

The device is intended for use in an electromagnetic environment in which radiated RF interference is controlled. The customer or user of the equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the equipment, as recommended below, depending on the maximum transmitting power of the communications equipment.

Power output maximum assigned value of the sender (W)	Separation distance depending on the frequency of the transmitter (m)		
	150kHz - 80MHz	80MHz - 800MHz	800MHz - 2.5GHz
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.737
1	1.167	1.167	2.330
10	3.690	3.690	7.368
100	11.67	11.67	23.300

For transmitters whose maximum rated transmit power is not given above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum transmit power characteristic of the transmitter in watts (W), according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the highest frequency range applies.  
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflections from structures, objects and people.








**ELECTRONIQUE DU MAZET**

**ZA ROUTE DE TENCE  
43520 LE MAZET SAINT VOY**

**Tél : +33 4 71 65 02 16  
Mail : [sav@electroniquedumazet.com](mailto:sav@electroniquedumazet.com)**

Your dealer / distributor :

A large, empty rounded rectangle with a dark blue border, intended for the user to write the name of their dealer or distributor.